

## Jailbreaker 1.0 Column Holder Adapter

### PST-ES

- Use any columns or spray emitters
- Plug-and-Play
- Comes with an HV cable for the user to connect to a liquid junction of their choice
- Magnetic coupler to insulated liquid junction fittings also available

#### Model number:

**PST-ES ; Basic adapter with HV cable**

**PST-ES-U: Adapter with an insulated uHPLC-compatible stainless steel union and magnetic coupler**

**PST-ES-Ru : Adapter with an insulated stainless steel reducing union and magnetic coupler. Compatible with a NanoViper® fitting from transfer tubing**

**PST-ES-T : Adapter with an insulated PEEK Tee and a gold wire and magnetic coupler**

**PST-ES-ST: Adapter with an insulated stainless steel Tee and magnetic coupler**

**PST-ES-SX: : Adapter with an insulated stainless steel cross and a gold wire and magnetic coupler.**

